


Anti-Pitrm1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q5Jrx3

Item Number: CM0006914



Introduction

PITRM1 Rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse PITRM1, a mitochondrial presequence protease (~117 kDa) involved in peptide degradation and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	PITRM1 Rabbit mAb
Remarks/Alias	MP1; PreP; SCAR30; PITRM1
Species	Human
GeneID (Human)	10531
GeneID	10531
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5JRX3
Protein Weight	117kDa
Shipping	Ice bag