


Anti-Pitrm1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Jrx3

Item Number: CM0030366



Introduction

Rabbit polyclonal antibody against human PITRM1, validated for WB, IF/ICC, and ELISA. Detects human, mouse, rat PITRM1. Ideal for mitochondrial peptide degradation and neurodegenerative disease research.

[Learn More](#)

Parameter	Details
Product Name	PITRM1 Rabbit pAb
Remarks/Alias	MP1; PreP; SCAR30; PITRM1
Species	Human
Gene ID (Human)	10531
Gene ID	10531
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 758-1037 of human PITRM1 (NP_055704.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5JRX3
Protein Weight	117kDa
Shipping	Ice bag