

Anti-Pituitary Homeobox 2 (Pitx2) Rabbit Monoclonal Antibody For Wb, Elisa - Q99697

Item Number: CM0006876



Introduction

Recombinant rabbit monoclonal antibody against human PITX2 transcription factor, validated for Western blot and ELISA applications. Ideal for developmental biology, left-right asymmetry, and disease mutation research.

[Learn More](#)

Parameter	Specification
Product Name	Pituitary homeobox 2 (PITX2) Rabbit mAb
Remarks/Alias	RS, RGS, ARP1, Brx1, IDG2, IGDS, IHG2, PTX2, RIEG, ASGD4, IGDS2, IRID2, Otlx2, RIEG1, Pituitary homeobox 2 (PITX2)
Species	Human
GeneID (Human)	5308
GeneID	5308
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Pituitary homeobox 2 (PITX2) (NP_700475.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99697
Protein Weight	35kDa
Shipping	Ice bag