

Anti-Pitx2 (Pituitary Homeobox 2) Rabbit Polyclonal Antibody For Wb, Elisa - Q99697

Item Number: CM0018941



Introduction

Rabbit polyclonal antibody against human PITX2 (Q99697), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying myogenesis, left-right asymmetry, and Peters anomaly.

[Learn More](#)

Parameter	Value
Product Name	Pituitary homeobox 2 (PITX2) Rabbit pAb
Remarks/Alias	RS; RGS; ARP1; Brx1; IDG2; IGDS; IHG2; PTX2; RIEG; ASGD4; IGDS2; IRID2; Otlx2; RIEG1; Pituitary homeobox 2 (PITX2)
Species	Human
GeneID	5308
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PITX2 (NP_000316.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99697
Protein Weight	35kDa
Shipping	Ice bag