

Anti-Pitx3 Rabbit Monoclonal Antibody For Wb, Elisa - 075364

Item Number: CM0003167

Introduction

PITX3 Rabbit Monoclonal Antibody validated for Western blot and ELISA, reacting with human, mouse, and rat. Suitable for research on dopaminergic neuron development and lens differentiation.

[Learn More](#)

Parameter	Value
Product Name	PITX3 Rabbit mAb
Remarks/Alias	ASMD; ASOD; PTX3; ASGD1; CTPP4; CTRCT11; PITX3
Species	Human
GeneID (Human)	5309
GeneID	5309
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PITX3 (075364).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075364
Protein Weight	32kDa
Shipping	Ice bag