


Anti-Piwil1 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q96J94

Item Number: CM0011057



Introduction

Rabbit monoclonal antibody against human PIWIL1 (HIWI, MIWI). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for germ cell and transposon silencing research.

[Learn More](#)

Attribute	Details
Product Name	PIWIL1 Rabbit mAb
Remarks/Alias	HIWI; MIWI; PIWI; CT80.1; PIWIL1
Species	Human
GeneID (Human)	9271
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human PIWIL1 (Q96J94).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96J94
Protein Weight	99kDa
Shipping	Ice bag