


Anti-Piwil1 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Jmb7

Item Number: CM0019447



Introduction

Rabbit polyclonal antibody against PIWIL1, validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for piRNA pathway and germ cell research.

[Learn More](#)

Parameter	Value
Product Name	PIWIL1 Rabbit pAb
Remarks/Alias	MIWI; PIWIL1
Species	Mouse
GeneID (Human)	9271
GeneID	57749
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of mouse PIWIL1 (NP_067286.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9JMB7
Protein Weight	99kDa
Shipping	Ice bag