

Anti-Piwil2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tc59

Item Number: CM0025190



Introduction

Rabbit anti-PIWIL2 polyclonal antibody validated for Western blot and ELISA, cross-reacts with mouse. Ideal for piRNA pathway, spermatogenesis, and germline integrity research. UniProt Q8TC59.

[Learn More](#)

Parameter	Value
Product Name	PIWIL2 Rabbit pAb
Remarks/Alias	CT80; HILI; mili; PIWIL1L; PIWIL2
Species	Human
GeneID (Human)	55124
GeneID	55124
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PIWIL2 (Q8TC59).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8TC59
Protein Weight	110kDa
Shipping	Ice bag