

Anti-Pjvk Rabbit Monoclonal Antibody For Wb, Elisa - Q0Zlh3

Item Number: CM0007220

Introduction

Anti-PJVK rabbit monoclonal antibody validated for WB and ELISA. Targets human pejvakin (DFNB59), a peroxisome-associated protein essential for hearing protection against noise-induced damage.

[Learn More](#)

Parameter	Value
Product Name	PJVK Rabbit mAb
Remarks/Alias	DFNB59; PJVK
Species	Human
GeneID	494513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-352 of human PJVK (NP_001036167.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q0ZLH3
Protein Weight	40kDa
Shipping	Ice bag