

Anti-Pjvk Rabbit Monoclonal Antibody For Wb, Elisa - Q0Zlh3

Item Number: CM0007198

Introduction

Anti-PJVK Rabbit Monoclonal Antibody validated for WB and ELISA. Recognizes human pejvakin (40 kDa), a peroxisome protein that protects auditory hair cells from noise-induced damage. Suitable for hearing loss research.

[Learn More](#)

Parameter	Value
Product Name	PJVK Rabbit mAb
Remarks/Alias	DFNB59; PJVK
Species	Human
GeneID (human)	494513
GeneID	494513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-352 of human PJVK (NP_001036167.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q0ZLH3
Protein Weight	40kDa
Shipping	Ice bag