


Anti-Pka C-Alpha (Prkaca) Monoclonal Antibody (Kd Validated) For Wb, If/Icc, Elisa - P17612

Item Number: CM0002385



Introduction

KD validated mouse monoclonal antibody targeting PKA C-alpha (PRKACA). Reacts with human, mouse, rat. Ideal for WB, IF/ICC, ELISA and cAMP-dependent kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PKA C-alpha (PRKACA) Mouse mAb
Remarks/Alias	CAFD1; PKACA; PPNAD4; A)
Species	Human
Gene ID (Human)	5566
Gene ID	5566
Immunogen	Recombinant protein of human [KD Validated] PKA C-alpha (PRKACA)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17612
Protein Weight	41kDa
Shipping	Ice bag