

Anti-Pka C-Beta (Prkacb) Rabbit Monoclonal Antibody For Wb - P22694

Item Number: CM0007616

Introduction

High-quality rabbit monoclonal antibody targeting human PKA C-beta (PRKACB). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying cAMP signaling and kinase function.

[Learn More](#)

Parameter	Value
Product Name	PKA C-beta (PRKACB) Rabbit mAb
Remarks/Alias	CAFD2; PKA C-beta (PRKACB)
Species	Human
GeneID (Human)	5567
Immunogen	Synthetic peptide corresponding to amino acids 251-351 of human PKA C-beta (PRKACB)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22694
Protein Weight	41kDa
Shipping	Ice bag