

Anti-Pka C-Beta (Prkacb) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22694

Item Number: CM0027759

Introduction

Rabbit polyclonal antibody targeting human cAMP-dependent protein kinase catalytic subunit beta (PRKACB) validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Useful for cAMP signaling and kinase research.

[Learn More](#)

Parameter	Value
Product Name	PKA C-beta (PRKACB) Rabbit pAb
Remarks/Alias	CAFD2; PKACB; PKA C-beta; PKA C-beta (PRKACB)
Species	Human
GenelD (Human)	5567
GenelD	5567
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-351 of human PKA C-beta (PRKACB) (NP_002722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22694
Protein Weight	41kDa
Shipping	Ice bag