

Anti-Pka Rii α (Prkar2A) Rabbit Monoclonal Antibody For Wb, Elisa - P13861

Item Number: CM0009398



Introduction

High-quality rabbit monoclonal antibody against human cAMP-dependent protein kinase type II- α regulatory subunit (PRKAR2A). Suitable for WB and ELISA, specifically reacts with human PRKAR2A. Ideal for cAMP signaling research.

[Learn More](#)

Parameter	Value
Product Name	PKA RII α (PRKAR2A) Rabbit mAb
Remarks/Alias	PKR2; PKA RII α (PRKAR2A)
Species	Human
GeneID	5576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 312-404 of human PKA RII α (PRKAR2A) (P13861)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13861
Protein Weight	46kDa
Shipping	Ice bag