

Anti-Pkc Alpha + Beta + Gamma Monoclonal Antibody For Wb, Elisa - P05771 / P17252 / P05129 / Q02156

Item Number: CM0006327



Introduction

Rabbit monoclonal antibody targeting PKC alpha, beta, and gamma. Suitable for WB and ELISA. Reacts with human, mouse, and rat species. Ideal for studying protein kinase C signaling, with a conserved immunogen design ensuring broad isoform recognition.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha + beta + gamma Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; PKC alpha + beta + gamma
Species	Human
GeneID (Human)	5578/5579/5582/5581
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-500 of human PKC alpha + beta + gamma (NP_002728.2/NP_002729.2/NP_002730.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05771/P17252/P05129/Q02156
Protein Weight	78kDa
Shipping	Ice bag