

Anti-Pkc Alpha Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P17252

Item Number: CM0002485



Introduction

Rabbit monoclonal antibody targeting human PKC alpha (P17252). Validated in WB, IHC-P, IP, ELISA with cross-reactivity to mouse and rat. Ideal for studying cell signaling, cancer, and cardiac biology.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; PKC alpha
Species	Human
GeneID (Human)	5578
GeneID	5578
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 573-672 of human PKC alpha (P17252).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag