


Anti-Pkc Alpha Rabbit Monoclonal Antibody For Wb/Elisa - P17252

Item Number: CM0007148



Introduction

High-quality rabbit monoclonal antibody targeting human PKC alpha. Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for cancer, cardiac hypertrophy, and cell signaling research. ~77kDa.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; PKC alpha
Species	Human
GeneID (Human)	5578
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PKC alpha. (NP_002728.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag