


Anti-Pkc Alpha Rabbit Polyclonal Antibody For Wb, Elisa - P17252

Item Number: CM0018711



Introduction

Rabbit polyclonal antibody targeting PKC alpha for WB and ELISA. Validated in human, mouse, rat. Recognizes a 77kDa protein. Ideal for studying kinase signaling in cancer, cardiovascular, and immune research.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit pAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; PKC alpha
Species	Human
Gene ID (Human)	5578
Gene ID	5578
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 573-672 of human PKC alpha (NP_002728.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag