


Anti-Pkc Alpha Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17252

Item Number: CM0014016



Introduction

High-quality rabbit polyclonal antibody against human PKC alpha, validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for research into cell signaling, cancer, and cardiovascular biology.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit pAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; PKC alpha
Species	Human
GeneID (Human)	5578
GeneID	5578
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 573-672 of human PKC alpha (NP_002728.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag