


Anti-Pkc Alpha Rabbit Polymab® For Wb And Elisa - P17252

Item Number: CM0007064



Introduction

High-specificity monoclonal antibody targeting human PKC alpha, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for cell signaling, cancer research, and kinase studies.

[Learn More](#)

Parameter	Value
Product Name	PKC alpha Rabbit PolymAb®
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI +/-; PKCalpha; PKC-alpha
Species	Human
GeneID (Human)	5578
GeneID	5578
Immunogen	Recombinant protein of human PKC alpha.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252
Protein Weight	77 kDa
Shipping	Ice bag