

Anti-Pkc-Beta Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05771

Item Number: CM0006653



Introduction

Offered by CamelBio, this rabbit monoclonal antibody against PKC-beta (PRKCB) is validated for WB, IHC-P, and ELISA. It cross-reacts with human, mouse, and rat, making it ideal for studying B-cell activation, apoptosis, and insulin signaling.

[Learn More](#)

Parameter	Value
Product Name	PKC-beta Rabbit mAb
Remarks/Alias	PKCB; PRKCB1; PRKCB2; PKCI(2); PKCbeta; PKC-beta
Species	Human
GeneID (Human)	5579
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 260-340 of human PKC-beta (NP_997700.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05771
Protein Weight	77kDa
Shipping	Ice bag