

Anti-Pkc-Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05771

Item Number: CM0015583



Introduction

Rabbit polyclonal antibody against human PKC-beta (PRKCB). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for signal transduction, cancer and immunology research.

[Learn More](#)

Parameter	Specification
Product Name	PKC-beta Rabbit pAb
Remarks/Alias	PKCB; PRKCB1; PRKCB2; PKCI(2); PKCbeta; PKC-beta
Species	Human
GeneID (Human)	5579
GeneID	5579
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 260-340 of human PKC-beta (NP_997700.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05771
Protein Weight	77kDa
Shipping	Ice bag