


Anti-Pkc Delta Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q05655

Item Number: CM0010260



Introduction

Rabbit monoclonal antibody targeting human PKC delta (Q05655). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PKC delta Rabbit mAb
Remarks/Alias	MAY1; PKCD; ALPS3; CVID9; nPKC-delta; PKC delta
Species	Human
GeneID (Human)	5580
GeneID	5580
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 547-653 of human PKC delta (Q05655)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05655
Protein Weight	78 kDa
Shipping	Ice bag