


Anti-Pkc Epsilon Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q02156

Item Number: CM0019960



Introduction

High-quality rabbit polyclonal antibody targeting PKC epsilon (PKCE/nPKC-epsilon). Validated for WB, IHC-P, ELISA and reactive across human, mouse, rat. Ideal for studying cell adhesion, migration, cancer invasion, and immune response.

[Learn More](#)

Parameter	Value
Product Name	PKC epsilon Rabbit pAb
Remarks/Alias	PKCE; nPKC-epsilon; PKC epsilon
Species	Human
GeneID (Human)	5581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-737 of human PKC epsilon (NP_005391.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02156
Protein Weight	84kDa
Shipping	Ice bag