


Anti-Pkc Gamma Monoclonal Antibody For Wb - P05129

Item Number: CM0005851



Introduction

High-specificity rabbit monoclonal antibody targeting human PKC gamma, validated for Western blot and ELISA. Cross-reacts with mouse and rat. A valuable tool for research into neuronal signaling, synaptic plasticity, pain, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	PKC gamma Rabbit mAb
Remarks/Alias	PKCC; PKCG; SCA14; PKCI(3); PKCgamma; PKC-gamma; PKC gamma
Species	Human
Gene ID (Human)	5582
Gene ID	5582
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 604-697 of human PKC gamma (P05129)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P05129
Protein Weight	78kDa
Shipping	Ice bag