

Anti-Pkc Gamma Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05129

Item Number: CM0027969



Introduction

Rabbit polyclonal antibody targeting PKC gamma, validated in WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for neuroscience and kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	PKC gamma Rabbit pAb
Remarks/Alias	PKCC; PKCG; SCA14; PKCI(3); PKCgamma; PKC-gamma; PKC gamma
Species	Human
GeneID (human)	5582
GeneID	5582
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 558-697 of human PKC gamma (NP_002730.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05129
Protein Weight	78kDa
Shipping	Ice bag