

Anti-Pkc Mu (Prkd1) Rabbit Monoclonal Antibody For Wb And Elisa - Q15139

Item Number: CM0007092

Introduction

High-specificity rabbit monoclonal antibody targeting human PKC mu (PRKD1, Q15139). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for kinase signaling, cell migration, Golgi transport, and cardiac research.

[Learn More](#)

Parameter	Value
Product Name	PKC mu Rabbit mAb
Remarks/Alias	PKD; PKCM; CHDED; PRKCM; PKC-MU; PKC mu
Species	Human
GeneID (Human)	5587
GeneID	5587
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PKC mu (Q15139).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15139
Protein Weight	102kDa
Shipping	Ice bag