

## Anti-Pkc Mu Rabbit Polyclonal Antibody For Wb And Elisa - Q15139

Item Number: CM0021798

### Introduction

Rabbit polyclonal antibody targeting human PKC mu (PRKD1) suitable for Western blot and ELISA. Cross-reacts with mouse. Useful for studies on cell signaling, cancer, and cardiac hypertrophy.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | PKC mu Rabbit pAb                                                                                      |
| Remarks/Alias    | PKD; PKCM; CHDED; PRKCM; PKC-MU; PKC mu                                                                |
| Species          | Human                                                                                                  |
| Gene ID (Human)  | 5587                                                                                                   |
| Gene ID          | 5587                                                                                                   |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 800-900 of human PKC mu (NP_002733.2) |
| Source           | Rabbit                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                  |
| Application      | WB, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse                                                                                           |
| SWISS            | Q15139                                                                                                 |
| Protein Weight   | 102kDa                                                                                                 |
| Shipping         | Ice bag                                                                                                |