

Anti-Pkc Zeta Monoclonal Antibody For Wb, Elisa - Q05513

Item Number: CM0007061

Introduction

CamelBio offers validated Anti-PKC zeta Rabbit mAb for WB, ELISA. Target protein kinase C zeta (PRKCZ, PKC2) functions in PI3K pathway, NF-kB activation, and cell polarity. Ideal for human research in signaling, inflammation, and cancer.

[Learn More](#)

Parameter	Value
Product Name	PKC zeta Rabbit mAb
Remarks/Alias	PRKCZ; PKC-ZETA; PKC2; protein kinase C zeta; PKC zeta
Species	Human
GeneID	5590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PKC zeta (NP_002735.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q05513
Protein Weight	46kDa/56kDa/67kDa
Shipping	Ice bag