

Anti-Pkc Zeta Rabbit Polyclonal Antibody For Wb, Elisa - Q05513

Item Number: CM0016651

Introduction

Rabbit polyclonal antibody against human PKC zeta (SWISS Q05513), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying PI3K, MAPK, and NF-kappa-B signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	PKC zeta Rabbit pAb
Remarks/Alias	PKC2; PKC-ZETA; PKC zeta
Species	Human
GeneID (Human)	5590
GeneID	5590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 476-547 of human PKC zeta (NP_002735.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05513
Protein Weight	68kDa
Shipping	Ice bag