


Anti-Pkce Rabbit Monoclonal Antibody For Wb, Elisa - Q02156

Item Number: CM0006392



Introduction

Rabbit monoclonal antibody targeting PKCε (nPKC-epsilon) for Western blot and ELISA. Shows reactivity across human, mouse, and rat species. Suitable for studies on kinase signaling in cell adhesion, migration, and immune response.

[Learn More](#)

Parameter	Value
Product Name	PKCε Rabbit mAb
Remarks/Alias	PKCE; nPKC-epsilon; PKCε
Species	Human
Gene ID	5581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-737 of human PKCε (NP_005391.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02156
Protein Weight	84 kDa
Shipping	Ice bag