

Anti-Pkd1L2 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z442

Item Number: CM0015835

Introduction

Rabbit polyclonal antibody against human PKD1L2, validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for studying polycystin-1-like protein 2 function.

[Learn More](#)

Parameter	Value
Product Name	PKD1L2 Rabbit pAb
Remarks/Alias	PC1L2; PKD1L2
Species	Human
GeneID (human)	114780
GeneID	114780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human PKD1L2 (NP_001265352.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q7Z442
Protein Weight	273kDa
Shipping	Ice bag