


Anti-Pkd2 (Polycystin-2) Rabbit Polyclonal Antibody For Wb, Elisa - Q13563

Item Number: CM0026971



Introduction

Rabbit polyclonal antibody targeting human PKD2 (Polycystin-2), suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for kidney disease and ion channel signaling research.

[Learn More](#)

Parameter	Value
Product Name	PKD2 Rabbit pAb
Remarks/Alias	PC2; PKD4; Pc-2; APKD2; TRPP2
Species	Human
Gene ID (human)	5311
Gene ID	5311
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 731-830 of human PKD2 (NP_000288.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13563
Protein Weight	110kDa
Shipping	Ice bag