


Anti-Pkd2 Rabbit Polyclonal Antibody For Wb, Elisa - Q13563

Item Number: CM0025753



Introduction

Rabbit polyclonal antibody targeting human Polycystin-2 (PKD2). Suitable for WB and ELISA; cross-reacts with human, mouse, rat. Ideal for studying calcium channel function and polycystic kidney disease.

[Learn More](#)

Parameter	Value
Product Name	PKD2 Rabbit pAb
Remarks/Alias	PC2; PKD4; Pc-2; APKD2; TRPP2; PKD2
Species	Human
GeneID (Human)	5311
GeneID	5311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 879-968 of human PKD2 (NP_000288.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13563
Protein Weight	110kDa
Shipping	Ice bag