


Anti-Pkm1/2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14618

Item Number: CM0018609



Introduction

Rabbit polyclonal antibody against PKM1/2 (SWISS P14618) validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for glycolysis and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PKM1/2 Rabbit pAb
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM1/2
Species	Human
Gene ID (Human)	5315
Gene ID	5315
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag