


Anti-Pkm1/2 Rabbit Polymab® Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P14618

Item Number: CM0008594



Introduction

High-quality rabbit monoclonal anti-PKM1/2 antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for glycolysis and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	PKM1/2 Rabbit PolymAb®
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30
Species	Human
GeneID (Human)	5315
GeneID	5315
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 321-420 of human PKM1/2 (NP_001193728.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag