


Anti-Pkm2 Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P14618

Item Number: CM0002995



Introduction

Rabbit monoclonal antibody specific to human PKM2. Validated for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat. Targets the key glycolytic enzyme and tumor biomarker pyruvate kinase M2. Ideal for cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PKM2-specific Rabbit mAb
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM2-specific
Species	Human
GeneID (Human)	5315
GeneID	5315
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 380-480 of human PKM2-specific (NP_002645.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag