


Anti-Pkm2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P14618

Item Number: CM0015422



Introduction

Rabbit polyclonal antibody against PKM2, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat PKM2. Suitable for glycolysis and cancer metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PKM2 Rabbit pAb
Remarks/Alias	PK3, TCB, p58, OIP3, PKM2, CTHBP, THBP1, HEL-S-30
Species	Human
GeneID	5315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-221 of PKM2 (NP_001193728.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag