


Anti-Pkm2-Specific Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P14618

Item Number: CM0021648



Introduction

Rabbit polyclonal antibody targeting human PKM2 (58kDa), validated for WB, IF/ICC, ELISA with cross-reactivity to mouse, rat. Suitable for studying PKM2's dual roles in glycolysis and transcription regulation in cancer.

[Learn More](#)

Parameter	Specification
Product Name	PKM2-specific Rabbit pAb
Remarks/Alias	PK3; TCB; p58; OIP3; PKM2; CTHBP; THBP1; HEL-S-30; PKM2-specific
Species	Human
GeneID (Human)	5315
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 380-480 of human PKM2-specific (NP_002645.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14618
Protein Weight	58kDa
Shipping	Ice bag