


## Anti-Pkmyt1 Polyclonal Antibody For Wb And Elisa - Q99640

Item Number: CM0020550



### Introduction

Anti-PKMYT1 Rabbit Polyclonal Antibody validated for western blot and ELISA, cross-reacting with human, mouse, and rat PKMYT1 kinase. An essential tool for cell cycle and mitotic regulation research.

[Learn More](#)

| Parameter        | Specification                                                                                                      |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PKMYT1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | MYT1; PPP1R126; PKMYT1                                                                                             |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 9088                                                                                                               |
| GeneID           | 9088                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 75-362 of human PKMYT1 (NP_004194.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q99640                                                                                                             |
| Protein Weight   | 55kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |