


Anti-Pkp2 Polyclonal Antibody For Wb - Q99959

Item Number: CM0027839



Introduction

Rabbit polyclonal antibody targeting PKP2 (Plakophilin-2). Validated for WB and ELISA in human, mouse, and rat. Immunogen aa 578-837. For desmosome biology, cardiac electrophysiology, and influenza A host restriction research.

[Learn More](#)

Parameter	Value
Product Name	PKP2 Rabbit pAb
Remarks/Alias	ARVD9; PKP2
Species	Human
Gene ID (Human)	5318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 578-837 of human PKP2 (NP_001005242.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99959
Protein Weight	97kDa
Shipping	Ice bag