

Anti-Pkr/Eif2Ak2 Monoclonal Antibody For Wb, If/Icc, Elisa - P19525

Item Number: CM0002728



Introduction

Anti-PKR/EIF2AK2 rabbit monoclonal antibody validated for WB, IF/ICC, ELISA in human samples. Detects interferon-induced dsRNA-activated kinase involved in antiviral innate immunity. Suitable for viral infection and stress response research.

[Learn More](#)

Parameter	Value
Product Name	PKR/EIF2AK2 Rabbit mAb
Remarks/Alias	PKR; PRKR; DYT33; LEUDEN; EIF2AK1; PPP1R83; PKR/EIF2AK2
Species	Human
GeneID (Human)	5610
GeneID	5610
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-551 of human PKR/EIF2AK2 (P19525).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P19525
Protein Weight	62 kDa
Shipping	Ice bag