

Anti-Pkr/Eif2Ak2 Rabbit Polyclonal Antibody For Wb, Elisa - P19525

Item Number: CM0024808



Introduction

Rabbit polyclonal antibody against human PKR/EIF2AK2, validated for WB and ELISA. Recognizes the interferon-induced dsRNA-activated kinase critical for antiviral innate immunity, stress response, and signal transduction. Suitable for research on viral infection, cancer, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	PKR/EIF2AK2 Rabbit pAb
Remarks/Alias	PKR; PRKR; DYT33; LEUDEN; EIF2AK1; PPP1R83; PKR/EIF2AK2
Species	Human
Gene ID (Human)	5610
Gene ID	5610
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human PKR/EIF2AK2 (NP_002750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19525
Protein Weight	62kDa
Shipping	Ice bag