


Anti-Pla2G2A Rabbit Monoclonal Antibody For Wb, Elisa - P14555

Item Number: CM0003397



Introduction

High-quality rabbit monoclonal antibody targeting human PLA2G2A (sPLA2-IIA), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for research on inflammation, host defense, and lipid signaling.

[Learn More](#)

Parameter	Value
Product Name	PLA2G2A Rabbit mAb
Remarks/Alias	MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2; PLA2G2A
Species	Human
GeneID	5320
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 65-144 of human PLA2G2A (P14555).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14555
Protein Weight	16 kDa
Shipping	Ice bag