


Anti-Pla2G2A Rabbit Polyclonal Antibody For Wb, Elisa - P14555

Item Number: CM0013138



Introduction

Rabbit polyclonal antibody against human PLA2G2A, a secreted phospholipase A2 involved in inflammation and antimicrobial defense. Validated for WB and ELISA; suitable for lipid metabolism and host defense research.

[Learn More](#)

Parameter	Value
Product Name	PLA2G2A Rabbit pAb
Remarks/Alias	MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2; PLA2G2A
Species	Human
GeneID (Human)	5320
GeneID	5320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-144 of human PLA2G2A (NP_001155201.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14555
Protein Weight	16kDa
Shipping	Ice bag