

Anti-Pla2G2D Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Unk4

Item Number: CM0026103

Introduction

Rabbit polyclonal antibody against human PLA2G2D (SPLASH), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying secretory phospholipase A2's anti-inflammatory and immunoregulatory roles.

[Learn More](#)

Parameter	Value
Product Name	PLA2G2D Rabbit pAb
Remarks/Alias	SPLASH; sPLA2S; PLA2IID; sPLA2-IID; PLA2G2D
Species	Human
GeneID (Human)	26279
GeneID	26279
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-145 of human PLA2G2D (NP_036532.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNK4
Protein Weight	17kDa
Shipping	Ice bag