

Anti-Pla2G2D Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Unk4

Item Number: CM0030011

Introduction

High-quality anti-PLA2G2D rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and rat PLA2G2D, a secretory phospholipase A2 involved in anti-inflammatory responses. Ideal for immunology and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	PLA2G2D Rabbit pAb
Remarks/Alias	SPLASH; sPLA2S; PLA2IID; sPLA2-IID; PLA2G2D
Species	Human
Gene ID (Human)	26279
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-145 of human PLA2G2D (NP_036532.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UNK4
Protein Weight	17kDa
Shipping	Ice bag