

Anti-Pla2G3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nz20

Item Number: CM0020616

Introduction

Rabbit polyclonal antibody targeting human PLA2G3 (sPLA2-III), validated for WB and ELISA. Reacts with human and mouse. PLA2G3 is a secretory phospholipase A2 involved in lipid metabolism, inflammation, and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	PLA2G3 Rabbit pAb
Remarks/Alias	SPLA2III; sPLA2-III; GIII-SPLA2; PLA2G3
Species	Human
GeneID (Human)	50487
GeneID	50487
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human PLA2G3 (NP_056530.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZ20
Protein Weight	57kDa
Shipping	Ice bag