


Anti-Pla2G4A Rabbit Polyclonal Antibody For Wb And Elisa - P47712

Item Number: CM0017327



Introduction

Rabbit polyclonal antibody against human cytosolic phospholipase A2 (PLA2G4A, cPLA2 α). Validated for WB and ELISA, detects human and mouse. Target involved in eicosanoid biosynthesis and inflammatory response. - P47712

[Learn More](#)

Parameter	Value
Product Name	PLA2G4A Rabbit pAb
Remarks/Alias	GURDP; cPLA2; PLA2G4; cPLA2-alpha; PLA2G4A
Species	Human
GeneID (Human)	5321
GeneID	5321
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 550-749 of human PLA2G4A (NP_077734.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47712
Protein Weight	85kDa
Shipping	Ice bag