

Anti-Pla2G4D/Cpla2Delta Polyclonal Antibody For If/Icc, Elisa - Q86Xp0

Item Number: CM0015742



Introduction

Rabbit polyclonal antibody targeting human cytosolic phospholipase A2 delta (PLA2G4D/cPLA2delta) for IF/ICC and ELISA. Involved in lipid metabolism and inflammation, expressed in stratified epithelia.

[Learn More](#)

Parameter	Value
Product Name	PLA2G4D Rabbit pAb
Remarks/Alias	cPLA2delta; PLA2G4D
Species	Human
Gene ID (Human)	283748
Gene ID	283748
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-680 of human PLA2G4D (NP_828848.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q86XP0
Protein Weight	92kDa
Shipping	Ice bag